

OCTOBER-DECEMBER 2024

# GOING GREEN



## Sustainability Campaign

Make a real impact by Going Green with NSK. Our Going Green campaign rescues used and non-performing handpieces and equipment from the landfill, giving them a second life through refurbishment and repair. Even items that can't be repaired are recycled responsibly. When you upgrade your equipment, you're not just enhancing your practice — you're contributing to a sustainable future.

Through our recycling and refurbishment programs, we support charitable organizations and uplift underserved communities both locally and globally. Join us in making a difference!



SCAN TO FIND OUT MORE



\*Selected products only.



1300 44 33 21 | [www.australia.nsk-dental.com](http://www.australia.nsk-dental.com) | [info@nskoceania.com.au](mailto:info@nskoceania.com.au)

# GO GREEN AND TRADE-IN YOUR HANDPIECES

Join us in our mission for a healthier planet and brighter smiles! Together, we can make a real difference. Go green for a sustainable future by trading in your old equipment and receiving up to 30% off on replacement items.\* We accept trade-ins of any make or model, even if they're not in working condition. However, disassembled handpieces, spare parts, or components cannot be accepted.

\*Discount applies to selected items only.



## Ti-Max Series



**MODEL Ti-Max Z900L**

ORDER CODE **P1111**

- Optic
- 26W power
- Standard head

**RRP \$1,981**

**TRADE-IN PRICE \$1,387**

**SAVE \$594**



**MODEL Ti-Max Z24L**

ORDER CODE **C1126001**

- 1:1 Direct Drive
- Single Spray
- Miniature head
- For CA burs (ø2.35)

**RRP \$1,491**

**TRADE-IN PRICE \$1,044**

**SAVE \$447**



**MODEL Ti-Max X65L**

ORDER CODE **H273**

- Clean Head System
- Single Spray
- For HP burs (ø2.35)

**RRP \$1,955**

**TRADE-IN PRICE \$1,369**

**SAVE \$587**



**MODEL Ti-Max X-DSG20L**

ORDER CODE **C1068**

- Optic
- 20:1 Reduction
- Dismantable

**RRP \$2,932**

**TRADE-IN PRICE \$2,052**

**SAVE \$880**

## Surgic Pro Optic



**MODEL Surgic Pro Optic Complete Set**

ORDER CODE **Y1001933**

- Control Unit
- SGL70M LED Micromotor FC-78 Foot Control
- X-SG20L Optic Handpiece (20:1 Reduction)
- Irrigation Tube
- 5 pcs and other accessories

**RRP \$8,726**

**TRADE-IN PRICE \$6,108**

**SAVE \$2,618**



## Varios-Combi Pro



**MODEL Varios-Combi Pro Basic Set**

ORDER CODE **Y1004194**

- Frequency: 28-32 kHz
- Dimensions: W 225 x D 290 x H 162
- Unit weight : 2.6 kg

**RRP \$7,856**

**TRADE-IN PRICE \$5,499**

**SAVE \$2,357**

**Get 2 FREE additional** Varios Combi Pro Prophy-Handpieces with the Varios Combi Pro Basic Set

## GOING GREEN CAMPAIGN



### FX Series Handpieces



**MODEL FX25M**  
ORDER CODE **C1052001**

- 1:1 Direct Drive
- For CA burs (ø2.35)
- Push Button Chuck
- Ball bearing
- Max speed: 40,000 min-1

**RRP \$462**

**TRADE-IN PRICE \$393**

**SAVE \$69**



**MODEL FX57M**  
ORDER CODE **C1054001**

- 4:1 Reduction
- For screw-in cups and brushes
- 360° Rotation
- Paste Defense System
- Max speed: 5,000 min-1

**RRP \$508**

**TRADE-IN PRICE \$432**

**SAVE \$76**



**MODEL FX65M**  
ORDER CODE **H1014001**

- For HP burs (ø2.35)
- Stainless Steel body
- Clean Head System
- Max speed: 40,000 min-1

**RRP \$461**

**TRADE-IN PRICE \$392**

**SAVE \$69**

### iCare

- Preset solution volume based on handpiece type (high-speed handpieces, low-speed handpieces)
- Chuck cleaning function
- Purging select mode
- Integrated large oil tank (1.2 Litres)



**MODEL iCare C2 Type**  
ORDER CODE **Y1002787**

- Dual Standard Ports + Dual Rotation Ports

**RRP \$5,489**

**TRADE-IN PRICE \$4,117**

**SAVE \$1,372**



**MODEL iCare C3 Type**  
ORDER CODE **Y1002788**

- Single Standard Ports + Triple Rotation Ports

**RRP \$5,863**

**TRADE-IN PRICE \$4,397**

**SAVE \$1,466**



### Volvere i7



**MODEL VOLVERE i7 RM**  
ORDER CODE **Y1003493**

- Dimensions: W 69 x D 185 x H 167 (mm)
- Weight: 930 g Standard Type Motor Handpiece
- Speed: 1,000-35,000 min-1

**RRP \$2,115**

**TRADE-IN PRICE \$1,586**

**SAVE \$529**

**MODEL VOLVERE i7 E**  
ORDER CODE **Y1003494**

- Dimensions: W 69 x D 185 x H 167 (mm)
- Weight: 930 g Standard Type Motor Handpiece
- Speed: 1,000-35,000 min-1

**RRP \$1,741**

**TRADE-IN PRICE \$1,306**

**SAVE \$435**

### ENDO-MATE DT2

**MODEL ENDO-MATE DT2 (230V)**  
ORDER CODE **Y1004230**

- Control Unit
- Motor and cord
- MP-F20R head (20:1)
- Handpiece stand
- AC adaptor

**RRP \$2,348**

**TRADE-IN PRICE \$1,761**

**SAVE \$587**



### Varios 970 Optic

**MODEL Varios 970 Optic Complete Set**  
ORDER CODE **Y1001168**

- Control Unit
- Handpiece and cord
- 3 Scaling Tips (G4,G6,G8)
- Foot Control
- 2 x Irrigation bottles
- 3 x Tip Wrench with torque limiter
- AC Power Cord
- Sterilisation Case

**RRP \$4,526**

**TRADE-IN PRICE \$3,395**

**SAVE \$1,132**



# iClave 24

ORDER CODE **92412B0100**

- External Dimensions: W 451 x H 492 x D 669 (mm)
- Chamber Dimensions: 265 x 475 (mm)
- Chamber Capacity: 24 L
- Net Weight: 55 kg
- Supply Voltage CE: 230 V - 50/60Hz
- Maximum Power Consumption: 2,350 W
- Cycle Water Consumption\*: 700 ~ 800 cc
- Tank Capacity (clean and used water): 5 litres each
- Sound Level: 60 dB(A)

\*Water consumption varies depending on device operating conditions (e.g. chamber temperature, cycle selected, etc.).

MODEL **iClave24 Autoclave with Printer**

RRP **\$13,970**

**TRADE-IN PRICE \$12,573**

**SAVE \$1,397**

MODEL **iClave24 Autoclave without Printer**

RRP **\$13,640**

**TRADE-IN PRICE \$12,276**

**SAVE \$1,364**



## Technology & Function



### Optimised load capacity

The iClave 24 offers a true 24-litre capacity with a 5-tray system, allowing up to 7 kg of load, maximising chamber space.



### Fast cycle, excellent performance

The high-efficiency steam generator with stainless steel coil ensures even heat distribution and rapid drying, enabling quick and effective Class B and Class S cycles.



### Internal ventilation

The internal ventilation system is durable, quiet, and ensures long-lasting, comfortable performance for the operator.



### A next-generation experience

The iClave 24's 4.3" capacitive color touchscreen is easy to use, even with gloves, making sterilisation management safe and efficient.



### A faster, safer, and more practical way to manage water

The iClave 24's chemical-resistant clean water tank is easy to fill, clean, and maintain, with quick connectors for safe and practical drainage.



### Easy traceability

The iClave 24 ensures safe sterilization cycles with multiple data tracking options, including Ethernet connection, USB transfer, and built-in or external printing for reports and labels.

All prices as listed. Valid until 31 December 2024. Not to be used in conjunction with any other NSK promotions. All prices are as quoted in AUD and include GST.

SCAN FOR  
FULL NSK  
ICLAVE 24  
RANGE

