



NSK

Outstanding Performance &
Excellent versatility



new


NSK presents **its latest**
Ultrasonic Scaler
Powered by iPiezo engine

Varios 970
iPiezo engine

Varios 970 iPiezo engine

NSK new technology, iPiezo engine ensures stable power output by automatically adjusting to the optimum vibration frequency depending on loading status, which offers you more effective treatment in all applications with wide range of tips.


Varios 970 well balanced slim and lightweight handpieces with LED



Competitor A

NSK

Slim ergonomically designed handpiece guarantees excellent accessibility, high visibility and comfortable operation.



LED Optic VA2-LUX-HP Optic Handpiece



Non-Optic VA2-HP Non-Optic Handpiece



NSK LED delivers natural daylight quality light bringing significant advantages to the way you work.

Clearer : LEDs generate natural daylight quality light which assures clearer vision when compared to halogen light.

Durable : LED light sources are longer lasting and more durable than halogen bulbs.

Safer : LED light is safer as it generates less heat than a halogen bulb. Natural daylight quality LED light is easy on the eye.

